

TabSafe

*Secure and simplified
medication management*

www.Healthcominc.com

TabSafe

The TabSafe system is the latest advancement for in-home medication management. Designed with a dual-locking feature, TabSafe offers enhanced security features and solutions for even the most complex dosing schedules. The system also delivers a flexible and easy to use programming module.



- Easy-to-use loading and dispensing
- Flexible scheduling
- Safety and security for users and caregivers
- Online reporting capabilities



Healthcom provides tools that help seniors live independently in their own home.

1600 West Jackson Street
Sullivan, Illinois 61951
800 525-6237 toll free phone
800 392-8957 toll free fax

Ruth & Janet

"Before TabSafe, I was always worrying that my mom wasn't taking her pills or taking the wrong dose. It was just so confusing with so many medications and directions for when to take them. Now, I have peace of mind that Mom will take exactly the right pills at the right time. She's happy to know she can live independently much longer."

Ruth



TabSafe is the easy-to-use, effective solution for managing medications yet can accommodate even the most complex schedules. The security features and online reporting mean peace of mind for the user, family and caregivers. Improving medication adherence can lead to increased wellness and independence.

Easy Installation, Loading and Daily Use


Plugs in to electrical power and standard phone lines. Each of the 4 cartridges can be loaded either with individual or multi-dose medications. Medications can be loaded into the system by the caregiver, pharmacist or user. TabSafe then transmits and downloads the pre-loaded medication schedule. When a medication is due, an alert reminds the user to press a button which then triggers the medication to dispense into a drawer for release.

Flexible

TabSafe can accommodate any schedule; regardless of the level of complexity. Based on the patient's medication regimen, each of the 4 cartridges contain either 16 or 32 easy to fill individual compartments. Specific cartridges can be allocated as "PRN Only" to allow the dispensing of medications on an "as needed" basis and thus eliminating the risk of overdose. Custom visual reminders can also be added for activities of daily living such as "Check Blood Pressure" or "Meal Time".

Live 24 Hr Operator Intervention

TabSafe can be programmed to call into the CareLink Support Center. This allows trained professionals to monitor medication compliance and to respond quickly if medications are not accessed at the correct time.



Safe and Secure

TabSafe is designed as a secure, locking unit. Each individual medication cartridge remains secure until the TabSafe unit is unlocked for reloading. The dispensing drawer also remains locked until the actual release of the medication occurs. Optional security features include the use of a PIN number to gain access to programming capabilities.

Reporting Features

TabSafe calls into the CareNet site daily to report: on time medication dispenses (including specific dosing times), missed doses, any refill alerts that are appropriate, early doses and PRN's taken (if any). The robust reporting capabilities provide valuable, life-saving insight to caregivers, family members and health care providers regarding medication compliance.

"At least 1.5 million Americans are injured every year by medication errors."

- Institute of Medicine panelist Albert Wu, a drug safety specialist at Johns Hopkins University